

Thinking

EFFECTO LATO
inserto
BACK



2 A Teschio e denti, 1489.

RL 19058 v (B. 41 r; K/P 42 v).

Il vacuo della cassa dell'occhio, e 'i vacuo dell'osso sostenitore della guancia, e quello del naso e della bocca sono d'eguale profondità e terminano sotto il senso comune [sede centrale delle sensazioni] per linea perpendicolare. E ciascuna d'esse vacuità ha tanto di profondità quant'è la terza parte del volto dell'omo, cioè dal mento a capegli.



BACK

A NEW SOLUTION
The Trunk Total Support

Accessories



Headrest Versa Capitis[®] HEAD SUPPORTS

MOUNTING SYSTEM



Mounting plate for Capitis Headrest with horizontal adjustments



CAPITIS CONFORT



CAPITIS ANATOMICO





CAPITIS CONFORT DIMENSIONS (expressed in centimeters)



ARTICLE	A	B	C	D	E	F
PMDC01PC	21	14	13	15	08	04
PMDC01MC	30	14	13	24	10	05
PMDC01GC	42	16	14	34	14	07

CAPITIS ANATOMICO DIMENSIONS (expressed in centimeters)




ARTICLE	A	B	C	D	E	F
PMDC02XPC	18	13	08	13	04	02
PMDC02PC	24	17	12	19	06	04
PMDC02MC	28	19	13	21	08	06
PMDC02GC	32	25	16	23	08	08

Harnesses Versa Fixatis[®] FASTENERS



DYNAMIC BUTTERFLY HARNESS



ARTICLE	A	B	C	D	E
PMD 16XS	11	18	17	05	29
PMD 16S	17	20	21	05	34
PMD 16M	20	20	26	05	46
PMD 16L	20	25	30	06	52



CHEST HARNESS WITH SHOULDERS RETRACTION PADS



ARTICLE	L	H
PMD 17S	14	16
PMD 17M	16	18
PMD 17L	23	23



TRUNK HARNESS WITH ELASTIC SHOULDER RETRACTORS AND CHEST BUCKLE CLOSURE

- H = Height of the pad
- L = Width of the pad

ARTICLE	H	L
PMD 21XS	25	3.8
PMD 21S	35	3.8
PMD 21M	38	5
PMD 21L	45	5



Why Effecto Lato ?


? ENCOURAGES TRUNK EXTENSION



**Crash tested
ISO 7176-19**



**3 modular
configurable inserts**



Multi-adjustable

CLINICAL INFORMATION

- Poor control of the trunk
- Pelvis-Trunk instability
- Severe deformities

TECHNICAL FEATURES

- 3 modular configurable inserts as standard
- Two layers of foam. Self-modeling polymeric foam and reticulated foam for ventilation
- Microclimatic ventilated cover
- Aluminum flat shell T shaped
- Mounting hardware kit with 4 attaching points

SPARE PARTS

- Cover
- Mounting hardware kit

ACCESSORIES

- Thoracic Supports
- Fixatis Range Harnesses
- Capitis Range Headrest
- Headrest Mounting Plate

TTS EMBRACING SHAPES





inserto EFFECTO LATO

3 modular configurable inserts

MODEL	Paediatric			Adult				
	3033	3338	3642	3945	4245	4250	4550	4855
Width (cm)	30	33	36	39	42	42	45	48
Height (cm)	33	38	42	45	45	50	50	55


• Maximum Load (referred to size 4855): **136 KG**

• Compatible with wheelchair tubing diameter from **19** up to **25 mm.** and from **¾"** up to **1"**.

Pushing lumbar contoured insert



Flat Shell T shaped



Lumbar support insert



Technical Information

Backrest Technical Information					Patients Measures	
SIZE (cm)	Shell Height (cm)	Mounting on wheelchairs (outside tubes) (cm)	Reclination (°)	Rotation (°)	Trunk width with thoracic supports (cm)	
30x33	33	from 27 up to 33	-7 up to +7	-5 up to +5	13-23	With B type thoracic support bracket
					18-28	With A type thoracic support bracket
33x38	38	from 30 up to 36			16-26	With B type thoracic support bracket
					21-31	With A type thoracic support bracket
36x42	42	from 33 up to 39			19-29	With B type thoracic support bracket
					24-34	With A type thoracic support bracket
39x45	45	from 36 up to 42			22-32	With B type thoracic support bracket
					27-37	With A type thoracic support bracket
42x45	45	from 39 up to 45			25-35	With B type thoracic support bracket
					30-40	With A type thoracic support bracket
42x50	50	from 39 up to 45	25-35	With B type thoracic support bracket		
			30-40	With A type thoracic support bracket		
45x50	50	from 42 up to 48	28-38	With B type thoracic support bracket		
			33-43	With A type thoracic support bracket		
48x55	55	from 45 up to 51	31-41	With B type thoracic support bracket		
			36-46	With A type thoracic support bracket		


Wheelchair tubes diameter: from 19 up to 25 mm. and from 3/4" up to 1"

Thoracic Supports

Accessories


PADS

FLAT PAD




SIZES (cm.)					
H	6	8	10	12	12
L	8	10	12	15	18

CONTOUR PAD




SIZES (cm.)					
H	6	8	10	12	
L	12	14	16	19	

FLAT PAD W/ARMPIT CUT



SIZES (cm.)					
H	8	10	12	12	
L	10	12	15	18	

CONTOUR PAD W/ARMPIT CUT



SIZES (cm.)					
H	8	10	12		
L	14	16	19		

HARDWARE





MANUFACTURED BY



Via Montagna Z.I. - Lotto 41
72023 Mesagne (BR) - ITALY
Tel. +(39) 0831 777840
Email: info@promedicare.it
www.promedicare.eu

CERTIFIED COMPANY WITH QUALITY
MANAGEMENT SYSTEM ISO 9001 / EN ISO 13485

CERTIFIED
ISO 9001
ISO 13485



DISTRIBUTED BY:



CE 2 YEAR WARRANTY

The manufacturer reserves the right to make changes or replacements without notice.
The contents, the form, the creation, the design, the image and the structure of this catalog are the property of Pro Medicare srl.
Any improper and unauthorised use is strictly forbidden. Any imitation and/or explicit or not explicit copy can be legally persecuted.
Cover works by artist Giuseppe Ciraci - www.giuseppeciraci.com - Any reproduction is explicitly forbidden.