

ORDER



Invoicing Details	
Company name:	_____
Delivery date required:	___ / ___ / ___
Delivery address:	_____
P. Code:	_____ Town: _____
Tel.:	_____ Fax: _____
Order #:	_____ Date of the order: ___ / ___ / ___

Shipping Details	
Ship to:	_____

CUSTOMER STAMP AND SIGNATURE	

Each model includes as standard: modular structural kit, pads, cover.

Pro Medicare Srl warrants the devices functionality for a maximum period of **24 months**, covering all manufacturing defects from the first commissioning and **12 months** on covers and wear parts. The warranty is valid if the device is used as indicated within the instruction manual.

Note: Any request for variation in dimensions, materials and type compared to the standard defined for the specific user does not have CE marking, therefore the professional user becomes the manufacturer and has the obligation to draw up the technical documentation of the device (this is valid for the sections of the sheet with the letter W).

Images are purely illustrative and may not fully reflect reality.



INSERTO Back CONCO (structural kit with padding and cover)															
Size (cm)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	W** <small>tailored</small>
	2415	2421	2815	2821	3215	3221	3227	3618	3624	3630	4018	4024	4030		
effective shell width	19	19	23	23	27	27	27	31	31	31	35	35	35	/	
effective shell height	15	21	15	21	15	21	27	18	24	30	18	24	30	/	
*effective wheelchair external tubes width	24-26	24-26	28-30	28-30	32-34	32-34	32-34	36-38	36-38	36-38	40-42	40-42	40-42	/	
n° of pieces															
PADDING: <input checked="" type="checkbox"/> polyurethane elastic foam <input type="checkbox"/> memory foam <input checked="" type="checkbox"/> standard <input type="checkbox"/> optional no extra charge															

*indicative values

**indicate effective shell width x effective shell height

INSERTO Back CONCO (structural kit with padding, cover and multi-tension strap system with H 5 cm single straps)															
Size (cm)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	W** <small>tailored</small>
	2415	2421	2815	2821	3215	3221	3227	3618	3624	3630	4018	4024	4030		
effective shell width	19	19	23	23	27	27	27	31	31	31	35	35	35	/	
effective shell height	15	21	15	21	15	21	27	18	24	30	18	24	30	/	
*effective wheelchair external tubes width	24-26	24-26	28-30	28-30	32-34	32-34	32-34	36-38	36-38	36-38	40-42	40-42	40-42	/	
n° of pieces															
H 5 cm single straps	3 straps	4 straps	3 straps	4 straps	3 straps	4 straps	5 straps	3 straps	4 straps	5 straps	3 straps	4 straps	5 straps		
PADDING: <input checked="" type="checkbox"/> polyurethane elastic foam <input type="checkbox"/> memory foam <input checked="" type="checkbox"/> standard <input type="checkbox"/> optional no extra charge															

*indicative values

**indicate effective shell width x effective shell height

INSERTO Back CONCO (structural kit with padding, cover and multi-tension strap system with the combination of a H 10 cm central double strap and H 5 cm single straps)

Size (cm)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	W** <i>tailored</i>
	2415	2421	2815	2821	3215	3221	3227	3618	3624	3630	4018	4024	4030		
effective shell width	19	19	23	23	27	27	27	31	31	31	35	35	35	/	
effective shell height	15	21	15	21	15	21	27	18	24	30	18	24	30	/	
*effective wheelchair external tubes width	24-26	24-26	28-30	28-30	32-34	32-34	32-34	36-38	36-38	36-38	40-42	40-42	40-42	/	
n° of pieces															
combination of H 10cm central double strap and H 5cm single straps	NO	YES	NO	YES	NO	YES	YES	NO	YES	YES	NO	YES	YES		

PADDING: polyurethane elastic foam memory foam *standard*
 optional no extra charge

*indicative values

**indicate effective shell width x effective shell height

SPARE PARTS - Cover

Size (cm)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	W** <i>tailored</i>
	2415	2421	2815	2821	3215	3221	3227	3618	3624	3630	4018	4024	4030		
n° of pieces															

SPARE PARTS - Padding

Size (cm)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	W** <i>tailored</i>
	2415	2421	2815	2821	3215	3221	3227	3618	3624	3630	4018	4024	4030		
n° of pieces															

SPARE PARTS - Multi-tension strap system with H 5 cm single straps

Size (cm)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	W** <i>tailored</i>
	2415	2421	2815	2821	3215	3221	3227	3618	3624	3630	4018	4024	4030		
n° of pieces															

SPARE PARTS - Multi-tension strap system with the combination of H 10 cm central double strap and H 5 cm single straps

Size (cm)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	W** <i>tailored</i>
	2415	2421	2815	2821	3215	3221	3227	3618	3624	3630	4018	4024	4030		
n° of pieces	/		/		/			/			/				

**indicate effective shell width x effective shell height