



Order  Quotation

Invoicing Details	Shipping Details
Company name: _____ Delivery date required: ___/___/___ Delivery address: _____ P. Code.: _____ Town: _____ Tel.: _____ Fax: _____ # Order: _____ Date of the Order: ___/___/___	Ship to: _____ _____ <p style="text-align: center;"><b>CUSTOMER STAMP AND SIGNATURE</b></p>

Pro Medicare srl grants **24 months** warranty on production defects from the date of the first put into service and **12 months** on covers and components subjected to wearing. Warranty is void if the product is used in any other way than stated in the Instruction Manual.

**Note:** Any request for variation in dimensions, materials and type compared to the standard does not have **CE** marking, therefore the professional user becomes the manufacturer and has the obligation to draw up the technical documentation of the device.

*Images are purely illustrative and may not fully reflect reality.*

Ermes Standard Configuration			
<ul style="list-style-type: none"> <li>• Seat adjustable in width, depth also asymmetrically in adduction and abduction</li> <li>• Kit of inserts for a customised positioning of the pelvis</li> <li>• Reclining and height adjustable backrest</li> <li>• Kit of inserts for a customised positioning of the trunk</li> <li>• Multi-adjustable and removable containment and positioning kit:                             <ul style="list-style-type: none"> <li>– pelvis (adjustable in width, depth, height, adduction/abduction)</li> <li>– arm support (adjustable in width, depth, height, angle, intra/extra rotation)</li> </ul> </li> </ul>			
Sizes of the user			
	1 RH	depth of the seat RH	_____
	1 LH	depth of the seat LH	_____
	2 RH	length of the leg RH	_____
	2 LH	length of the leg LH	_____
	3 RH	depth of the thorax RH	_____
	3 LH	depth of the thorax LH	_____
	4	height of the shoulders from the seat plane	_____
	5	width of the thorax	_____
6	width of the pelvis	_____	
Ermes			
<input type="radio"/> size 2	color	<input type="radio"/> blue/grey	<input type="radio"/> fuchsia/grey
<input type="radio"/> size 3	color	<input type="radio"/> blue/grey	<input type="radio"/> fuchsia/grey
<input type="radio"/> size 4	color	<input type="radio"/> wine/grey	
Note: <input type="radio"/> = optional without up-charge <input type="checkbox"/> = optional with up-charge			

## Type of the seat

- MIDI Seat**  
lateral/anterior control  
(standard configuration)



- DEEP Seat**  
total control (PTS)  
optional without  
up-charge



- Greater closure of the seat-backrest angle ( $\pm 5^\circ$ )

## Thoracic Supports

- Flat pads

	Size (cm)				
	H	1	2	3	4
	L	7	7	8	10
	L	7.5	8.5	10	12

- Contoured pads

	Size (cm)				
	H	1	2	3	4
	L	6	7	8	10
	L	9.5	11	13	16

Side	Size	Bracket type*	Bracket shape
<input type="checkbox"/> RH		<input type="checkbox"/> Fixed <input type="checkbox"/> Swing-Away	<input type="checkbox"/> A <input type="checkbox"/> B
<input type="checkbox"/> LH		<input type="checkbox"/> Fixed <input type="checkbox"/> Swing-Away	<input type="checkbox"/> A <input type="checkbox"/> B

Side	Size	Bracket type*	Bracket shape
<input type="checkbox"/> RH		<input type="checkbox"/> Fixed <input type="checkbox"/> Swing-Away	<input type="checkbox"/> A <input type="checkbox"/> B
<input type="checkbox"/> LH		<input type="checkbox"/> Fixed <input type="checkbox"/> Swing-Away	<input type="checkbox"/> A <input type="checkbox"/> B

\*In the technical information at the end of this form you can find the dimensions reachable with the two types of brackets (A type = L/90° Bracket; B type = 2.5 cm Offset Bracket)

## Belts and Harnesses (values expressed in cm, see details in the technical information)

- Chest pad

 (with thoracic supports only)	Size	HxL
	<input type="checkbox"/> 1	(8x10)
	<input type="checkbox"/> 2	(9x13)
	<input type="checkbox"/> 3	(10x13)

- Dynamic Abduction Harness

	<input type="checkbox"/>	XS (2)*
	<input type="checkbox"/>	S (3)*
	<input type="checkbox"/>	M (4)*

- Pelvic belt

	Size	HxL
	<input type="checkbox"/> 1	(6x11)
	<input type="checkbox"/> 2	(7x13)
	<input type="checkbox"/> 3	(7x15)

- Dynamic Butterfly Harness

	<input type="checkbox"/>	XS (2)*
	<input type="checkbox"/>	S (3)*
	<input type="checkbox"/>	M (4)*

- Jacket Harness

	<input type="checkbox"/>	XS (2)*
	<input type="checkbox"/>	S (3)*
	<input type="checkbox"/>	M (4)*

- Contoured Pelvic Belt

	Size	HxL
	<input type="checkbox"/> 1	(6x11)
	<input type="checkbox"/> 2	(8x15)
	Value H is referred to the height of the central part of the pad	

\*The number in the brackets is referred to the relevant Ermes size. Therefore it is **advisable** to combine the default size.

## Hardrest with multi-axial joint and multi-adjustable arm

- Confort** with  
 detachable hardware  
 fixed hardware

	<input type="checkbox"/>	S
	<input type="checkbox"/>	M
	<input type="checkbox"/>	G

- Anatomic** with  
 detachable hardware  
 fixed hardware

	<input type="checkbox"/>	XS
	<input type="checkbox"/>	S
	<input type="checkbox"/>	M

- 3 piece pads supportive headrest** with  
 detachable hardware  
 fixed hardware

	<input type="checkbox"/> 1	(7x13)
	<input type="checkbox"/> 2	(9x20)
	values expressed in cm of the occipital pad only	

## Other accessories



**Knee elevating legrest**  RH  LH  
(with cover for angle knee)

separate footplates:  RH  LH

**Footplate options:** slotted separate footplates:  RH  LH

footboard (it does not allow asymmetry, adduction and abduction)



Shoes Holders



**Calf support**  RH  LH

Calf panel (it does not allow asymmetry, adduction and abduction)

**Lateral supporting pad for legs**  RH  LH

Padded anterior/posterior Calf strap

Padding for separate footplates

Padding for footboard

**Elbow lateral-posterior support**  RH  LH

**Grey detachable Tray**  
with body contouring recess and edges


standard version

larger version for daily activities

**Transparent detachable Tray**  
with body contouring recess and edges  
(only standard version)

**N.B.:** the dimensions of the tray in the standard grey color version and those of the transparent tray are the same, the naming changes (the size 2 of the grey tray will be the size 1 of the transparent tray, 3->2, 4->3)

standard version dimensions



**Size (cm)**

**L**

**P**


**P1**

**2** (for Ermes size 2) 40 35 12

**3** (for Ermes size 3) 48 40 13.5

**4** (for Ermes size 4) 58 45 15.5

larger version for daily activities dimensions



**Size (cm)**

**L**

**P**

**P1**

**1** (for Ermes size 2) 50 45 13.7

**2** (for Ermes size 3) 60 50 15.4

**3** (for Ermes size 4) 68 55.6 19.3

**Padding of the tray**  upper  lower

## Interface kits for the frames

(compatibility with other bases not manufactured by Pro Medicare Srl is under the responsibility of the professional user)



**Male interface plate** (to mount under the Ermes seat)



**Female interface plate** (to mount on other frames)



**Solid base for interface plate fixing**

In this case, please indicate: Width \_\_\_\_\_ ; Depth \_\_\_\_\_ ;

NOTES for the positioning \_\_\_\_\_

**Note:** the base does not include the kit for mounting on the frame; for special requests please contact the Technical-Sales Department



**Complete male/female interface plate**

