


**ORDER FORM**

<u>Invoicing Details</u>	<u>Shipping Details</u>
Company name: _____	Ship to: _____ _____ <b>CUSTOMER STAMP AND SIGNATURE</b>
Delivery date required: ___ / ___ / ___	
Delivery address: _____	
P. Code.: _____ Town: _____	
Tel.: _____ Fax: _____	
# Order: _____ Date of the Order: ___ / ___ / ___	

Pro Medicare srl grants **24 months** warranty on production defects from the date of the first put into service and **12 months** on covers and components subjected to wearing. Warranty is void if the product is used in any other way than stated in the Instruction Manual.

Note: Any request for variation in dimensions, materials and type compared to the standard does not have CE marking, therefore the professional user becomes the manufacturer and has the obligation to draw up the technical documentation of the device.

Images are purely illustrative and may not fully reflect reality.

<b>Standard Configuration</b>	
<ul style="list-style-type: none"> <li>• Frame Adacta Liga Light anodized satin colour Titanium</li> <li>• Fixed folding down back posts</li> <li>• Reclining mechanism by predetermined plate w/pin 3°, 9°, 15°, 22°</li> <li>• Continuous Tilt in Space by drive lever from 4° up to 30° (from 6° up to 39°)<sup>1</sup> <small><sup>1</sup>Tilt in space can be up-graded of 4° with an increment of the starting position of the frame by operating on the positioning holes of the caster fork</small></li> <li>• Adjustment of the seat depth</li> <li>• Height adjustment: High (49 cm from the floor) /low (43 cm from the floor)*</li> <li>• Detachable footplates with 75° fixed knee angle</li> <li>• Fixed Anti-tip system</li> <li>• Wheels Kit with brake</li> <li>• Rigid wooden seat base covered leather flex</li> <li>• Kit of ABS laterals with arm supports</li> <li>• Versa Inserto Trunk/Pelvis positioning system</li> <li>• Capitis Headrest</li> <li>• Fixatis Pelvic Belt</li> </ul>	 <p>The images are purely indicative and may not reflect the standard configuration</p>



◊ = Standard    ○ = Optional without up-charge    ◻ = Optional with up-charge  
If a selection is not made, the product will be delivered in its standard configuration

Frame sizes (Outside tubes Width x Depth cm)						Frame Type
<input type="radio"/> XS1 (38x46)	<input type="radio"/> S (41x46)	<input type="radio"/> M (41x50)	<input type="radio"/> M1 (44x50)	<input type="radio"/> L (47x46)	<input type="radio"/> XL (47x50)	◊ Rigid <input type="radio"/> Folding
<input type="checkbox"/> Pre-assembly of postural components (collars for flat base, backrest attachments, headrest hardware)						
<b>Knee Angle</b>	◊ 75° <input type="radio"/> 90° <input type="checkbox"/> vertical elevating and telescopic					
<b>Castors</b>	◊ 150mm (6") poly grey <input type="radio"/> black <input type="radio"/> 175mm (7") poly <input type="radio"/> grey <input type="radio"/> black					
<b>Posterior Wheels</b> (height 43 cm from seat to the floor with 300 mm and 500 mm wheels; height 47 cm from seat to the floor with 600 mm wheels)						
◊ 300mm (12") poly grey with parking brake lever						
<input type="radio"/> 300mm (12") poly black with parking brake lever						
<input type="radio"/> 300mm (12") pneumatic grey with parking brake lever						
<input type="checkbox"/> 300mm (12") poly black with drum brake						
<input type="checkbox"/> 300mm (12") poly grey with drum brake						
<input type="checkbox"/> 300mm (12") pneumatic grey with drum brake						
<input type="checkbox"/> 300mm (12") poly black with drum brake and parking brake lever						
<input type="checkbox"/> 300mm (12") poly grey with drum brake and parking brake lever						
<input type="checkbox"/> 300mm (12") pneumatic grey with drum brake and parking brake lever						
<input type="checkbox"/> 400mm (16") poly black with parking brake lever (only 45 cm height)						
<input type="checkbox"/> 400mm (16") poly black with drum brake (only 45 cm height)						
<input type="checkbox"/> 400mm (16") poly black with drum brake and parking brake lever (only 45 cm height)						
<input type="checkbox"/> 500mm (20") poly grey without handrims with parking brake lever						<input type="checkbox"/> handrims for 500mm (20") wheels
<input type="checkbox"/> 500mm (20") pneumatic grey without handrims with parking brake lever						
<input type="checkbox"/> 500mm (20") pneumatic grey with handrims with drum brake						
<input type="checkbox"/> 500mm (20") pneumatic grey with handrims with drum brake and parking brake lever						
<input type="checkbox"/> 600mm (24") pneumatic grey with handrims with parking brake lever (NOTE. The tilt in space is reduced to 20°, only 47 cm height)						



<b>Anti-tip system</b>	<input checked="" type="checkbox"/> <i>fixed</i> <input type="checkbox"/> detachable	
<b>Tilt in space</b>	<input checked="" type="checkbox"/> <i>with lever</i> <input type="checkbox"/> with pedal	
<b>Pelvis Positioning System</b>		
<input checked="" type="checkbox"/> <i>Inserto Seat Basico</i>		<input type="checkbox"/> Inserto Seat Mediano
<input type="checkbox"/> Inserto Seat Domino		
<input type="checkbox"/> Waterproof seat cover with elastic band		
<b>Trunk Positioning System</b>		
<input checked="" type="checkbox"/> <i>Inserto Back Effecto Lato</i>	<input checked="" type="checkbox"/> Lumbar pushing insert <input checked="" type="checkbox"/> RH <input checked="" type="checkbox"/> LH <input type="checkbox"/> Lumbar containment insert <input type="checkbox"/> RH <input type="checkbox"/> LH	
<input type="checkbox"/> Inserto Back Supporto		
<b>Backpost (Note: Push handles height 63 cm)</b>		
<input checked="" type="checkbox"/> <i>Kit for fixed push handles</i> <input type="checkbox"/> Kit for push handles adjustable in height by pin (max adjustment 10 cm with 2.5 cm gap) <input type="checkbox"/> Kit for telescopic push handles with lever block (max adjustment 10 cm) <input type="checkbox"/> Kit for angle adjustable single handlebar for rigid frame <input type="checkbox"/> Kit for telescopic and angle adjustable single handlebar <input type="checkbox"/> Quick-release pins for easier backrest inclination ( <i>not valid for crash test</i> )		








## Accessories


**Thoracic Supports for Inserto Back Supporto** *Please, indicate the type, the shapes and the size desired*

Type of Bracket	Shape of the Pad					
Fixed <input type="checkbox"/> RH <input type="checkbox"/> LH	Flat <input type="checkbox"/> RH <input type="checkbox"/> LH 	cm	<input type="checkbox"/> size 2	<input type="checkbox"/> size 3	<input type="checkbox"/> size 4	<input type="checkbox"/> size 5
		H	8	10	12	12
		L	8.5	10.5	13.5	16.5
Swing-Away <input type="checkbox"/> RH <input type="checkbox"/> LH	Contoured <input type="checkbox"/> RH <input type="checkbox"/> LH 	cm	<input type="checkbox"/> size 2	<input type="checkbox"/> size 3	<input type="checkbox"/> size 4	
		H	8	10	12	
		L	12.5	14.5	17.5	

**Thoracic Supports for Inserto Back Effecto Lato** *Please, indicate the type, the shapes and the size desired*

Type of Bracket	Shape of the Pad					
Fixed <input type="checkbox"/> RH <input type="checkbox"/> LH	Flat <input type="checkbox"/> RH <input type="checkbox"/> LH 	cm	<input type="checkbox"/> size 2	<input type="checkbox"/> size 3	<input type="checkbox"/> size 4	<input type="checkbox"/> size 5
		H	8	10	12	12
		L	10	12	15	18
Swing-Away <input type="checkbox"/> RH <input type="checkbox"/> LH	Contoured <input type="checkbox"/> RH <input type="checkbox"/> LH 	cm	<input type="checkbox"/> size 2	<input type="checkbox"/> size 3	<input type="checkbox"/> size 4	
		H	8	10	12	
		L	14	16	19	
Detachable <input type="checkbox"/> RH <input type="checkbox"/> LH	<b>Bracket shape</b> * <i>(for the selection see Technical Information)</i> <b>type A:</b> <input type="checkbox"/> RH <input type="checkbox"/> LH <b>type B:</b> <input type="checkbox"/> RH <input type="checkbox"/> LH					

<b>Pedal for barriers</b> <input type="checkbox"/> RH <input type="checkbox"/> LH  <input type="checkbox"/> <b>Vehicle Transportation Kit*</b> (Only frame height 43cm and with 300mm and 500mm posterior wheels) <i>*according to the ISO 7176-19 standard: do not use the wheelchair as a seat in the vehicle if the weight of the user is lower than 22 kg.</i>	<b>Padding for ABS Laterals</b> <input type="checkbox"/> RH <input type="checkbox"/> LH  <b>Abductor wedge</b> <i>only with rigid seat base</i>  <table border="1" style="display: inline-table; vertical-align: top;"> <tr> <td>cm</td> <td><input type="checkbox"/> S</td> <td><input type="checkbox"/> M</td> <td rowspan="3"> <input type="checkbox"/> removable with lever  </td> </tr> <tr> <td>H</td> <td>7.5</td> <td>9</td> </tr> <tr> <td>D</td> <td>9</td> <td>10</td> </tr> <tr> <td>L</td> <td>7.5</td> <td>9</td> <td> <input type="checkbox"/> removable with mounting block  </td> </tr> </table> <p><b>Note:</b> if not indicate the abductor wedge will be removable with lever</p> <p><b>Please indicate the position of holes under the seat</b>          (if not indicated the abductor wedge will be fixed in the standard position)  <input type="checkbox"/> standard  <input type="checkbox"/> not centred    <input type="checkbox"/> on the RH    <input type="checkbox"/> on the LH</p>	cm	<input type="checkbox"/> S	<input type="checkbox"/> M	<input type="checkbox"/> removable with lever 	H	7.5	9	D	9	10	L	7.5	9	<input type="checkbox"/> removable with mounting block 
cm	<input type="checkbox"/> S	<input type="checkbox"/> M	<input type="checkbox"/> removable with lever 												
H	7.5	9													
D	9	10													
L	7.5	9	<input type="checkbox"/> removable with mounting block 												
<b>Lower limb accessories</b> <input type="checkbox"/> <b>Separate Calf Supports swing-away</b> <input type="checkbox"/> <b>Padded Calf Strap</b> <b>Cover for legrest tube</b> <input type="checkbox"/> RH <input type="checkbox"/> LH <b>Cover for Knee Angle joint</b> <input type="checkbox"/> RH <input type="checkbox"/> LH <input type="checkbox"/> <b>Padding and cover of footplate in leatherette</b>	<b>Tray</b> <input type="checkbox"/> <b>Rectangular Wooden Tray w/lateral flip hardware</b> <input type="checkbox"/> <b>Rectangular Wooden Tray w/detachable hardware</b> <input type="checkbox"/> <b>Padding for Tray</b>														
<b>Capitis Headrest Comfort size:</b> <input type="checkbox"/> Small <input checked="" type="checkbox"/> <i>Medium</i> <input type="checkbox"/> Large															
<b>Fixatis Pelvic belt 2P.2A. with squeeze plastic buckle size:</b> <input type="checkbox"/> S <input checked="" type="checkbox"/> <i>M</i> <input type="checkbox"/> L															
<b>Fixatis Dynamic Butterfly Harness</b> <input type="checkbox"/> S <input type="checkbox"/> M <input type="checkbox"/> L															

Sizes armrest and lateral containment (cm)				
	Armrest		ABS Lateral	
	H	5	H	13
	L	33	L	28

## Technical Features

**(Note: The length dimensions, referred to the average values, are subjected to variations due to the displacement of the wheel plates. The values reported are intended for the standard configuration)**

Zeus Size	Maximum Load (kg)	Frame Weight (kg)	Heaviest component weight (Supporto Back) (kg)	Maximum overall dimensions in width/length (cm)						Pivot Width (cm)	Max Slope	Seat rail to floor (cm)		
				standard 300mm wheels	300mm wheels w/drum brake	standard 400mm wheels	400mm wheels w/drum brake	standard 500mm wheels	500mm wheels w/drum brake				600mm wheels	
<b>XS1</b>	<i>Rigido</i>	75	19.2	3.6	56/100	57/100	58/100	61/100	58/106	61/106	58/112	75	7°	43/49 with 300/500mm posterior wheels
	<i>Chiudibile</i>	75	19.5											
<b>S</b>	<i>Rigido</i>	75	19.5	3.9	59/100	60/100	61/102	64/102	61/106	64/106	61/112	75		
	<i>Chiudibile</i>	75	19.9											
<b>M</b>	<i>Rigido</i>	100	19.8	3.9	59/103	60/103	61/104	64/104	61/109	64/109	61/116	78		
	<i>Chiudibile</i>	75	20.2											
<b>M1</b>	<i>Rigido</i>	100	19.9	4.4	62/103	63/103	64/104	67/104	64/109	67/109	64/116	78		
	<i>Chiudibile</i>	75	20.3											
<b>L</b>	<i>Rigido</i>	100	20.1	5.1	65/100	66/100	67/102	70/102	67/106	70/106	67/112	75		
	<i>Chiudibile</i>	75	20.5											
<b>XL</b>	<i>Rigido</i>	100	20.3	5.1	65/103	66/103	67/104	70/104	67/109	70/109	67/116	78		
	<i>Chiudibile</i>	75	20.8											

Zeus Size	User Parameters		Frame Depth Adjustments (cm)	Cushion Size (W x D) (cm)	Supporto Back sizes (cm)				Effetto Lato Back sizes (cm)			
	Pelvis Width (cm)	Seat Depth (cm)			Backrest	Height	Shell Width	Trunk Width	Backrest	Height	Shell Width	*Trunk Width
<b>XS1</b>	Adjustable with inner inserts	35-45	from 44 to 52	40x45	3944	44	36	26-34	3945	45	39	w/bracket A 27-38 w/bracket B 22-32
<b>S</b>		35-45	from 44 to 52	42x45	4248	48	39	29-37	4245	45	42	w/bracket A 30-40 w/bracket B 25-35
<b>M</b>		40-50	from 48 to 56	42x50	4248	48	39	29-37	4250	50	42	w/bracket A 30-40 w/bracket B 25-35
<b>M1</b>		40-50	from 48 to 56	45x50	4550	50	42	32-40	4550	50	45	w/bracket A 33-43 w/bracket B 28-38
<b>L</b>		35-45	from 44 to 52	48x45	4852	52	45	35-43	4855	55	48	w/bracket A 36-46 w/bracket B 31-41
<b>XL</b>		40-50	from 48 to 56	48x50	4852	52	45	35-43	4855	55	48	w/bracket A 36-46 w/bracket B 31-41