

### Eolo High chair Technical Sheet

<p><b>Manufacturer identification</b></p>	<p>Pro Medicare s.r.l. Via Montagna Z.I. Lotto 41 72023 Mesagne (Br) P.IVA 01803920741 tel. 00 39 0831 777840 Website: <a href="http://www.promedicare.eu">www.promedicare.eu</a> e-mail: info@promedicare.it</p>																																													
<p><b>Device commercial name</b></p>	<p><b>EOLO High chair</b></p>																																													
<p><b>UDI-DI BASIC</b></p>	<p><b>805571347ERMESTIPCOKM</b></p>																																													
<p><b>Device code</b></p>	<p><b>D04-X01-Y</b> where X represents the color; if X = 0 blue/green pastel, if X = 6 purple/dove grey; Y represents size 1, 2, 3</p>																																													
<p><b>Code of classification (CND)</b></p>	<p><b>Y1809 HIGH CHAIRS/POSITIONING SYSTEM/SITTING POSITIONING SYSTEMS</b></p>																																													
<p><b>Medical device class</b></p>	<p>The Device is CE marked according to EU Regulation 2017/745 and it is a non sterile Class I Medical Device</p>																																													
<p><b>Intended use</b></p>	<p>EOLO is the multifunctional tilting high chair; it can be used by user with limited mobility to maintain sitting position. It includes: - an adjustment mechanism to adapt the height of the seat as needed; - a seat tilting system; - a seat/backrest angle reclining system.</p>																																													
<p><b>Technical features</b></p>	<p><b>STANDARD CONFIGURATION</b></p> <ul style="list-style-type: none"> <li>• Frame for Positioning System with multi-adjustable push handle, 4 swivel castors (2 with locks); high-low height adjustable and tilting by gas spring (-5° to +30°)</li> <li>• Anatomical seat adjustable in width and depth</li> <li>• Anatomic backrest adjustable in height with extension, and reclining by gas spring (-3° to +20°)</li> <li>• Armrests adjustable in width, height, depth and tilt</li> <li>• Side guides adjustable in width, height, depth and tilt</li> <li>• Fixed Thoracic supports adjustable in height, depth and tilt</li> <li>• Tilting wooden tray with recess and edges</li> <li>• Footboard with containment adjustable in height, tilt, and flexion-extension</li> </ul> <p><b>DIMENSIONS/WEIGHT/LOAD:</b></p> <table border="1" data-bbox="491 1451 1441 1854"> <thead> <tr> <th colspan="2"></th> <th>Size 1</th> <th>Size 2</th> <th>Size 3</th> </tr> </thead> <tbody> <tr> <td rowspan="9" style="writing-mode: vertical-rl; transform: rotate(180deg);"><b>EOLO</b></td> <td>Seat width (cm)</td> <td>17-24</td> <td>22-29</td> <td>27-34</td> </tr> <tr> <td>Seat depth (cm)</td> <td>19-26</td> <td>22-34</td> <td>31-43</td> </tr> <tr> <td>Trunk height (cm)</td> <td>26-38</td> <td>33-45</td> <td>41-53</td> </tr> <tr> <td rowspan="2">Trunk width (cm)</td> <td>bracket A 19-25</td> <td>bracket A 22-28</td> <td>bracket A 27-33</td> </tr> <tr> <td>bracket B 14-20</td> <td>bracket B 17-23</td> <td>bracket B 22-28</td> </tr> <tr> <td>Distance Popliteal Cable-Footboard (cm)*</td> <td>10-21</td> <td>10-26</td> <td>10-36</td> </tr> <tr> <td>Armrest height (cm)</td> <td>14-19</td> <td>20-26</td> <td>20-26</td> </tr> <tr> <td>Reclining backrest (°)</td> <td>-3+20</td> <td>-3+20</td> <td>-3+20</td> </tr> <tr> <td>Eolo max load (kg)</td> <td>25</td> <td>35</td> <td>50</td> </tr> <tr> <td>Eolo weight (kg)</td> <td>18</td> <td>21</td> <td>26</td> </tr> </tbody> </table> <p><b>PERFORMANCE AND DURABILITY:</b> If used in safe conditions, performance and durability of the structural parts is 5 years, while for the seat parts (padding) is 2 years; it is valid under normal conditions of use and proper maintenance as specified in the instruction manual.</p>			Size 1	Size 2	Size 3	<b>EOLO</b>	Seat width (cm)	17-24	22-29	27-34	Seat depth (cm)	19-26	22-34	31-43	Trunk height (cm)	26-38	33-45	41-53	Trunk width (cm)	bracket A 19-25	bracket A 22-28	bracket A 27-33	bracket B 14-20	bracket B 17-23	bracket B 22-28	Distance Popliteal Cable-Footboard (cm)*	10-21	10-26	10-36	Armrest height (cm)	14-19	20-26	20-26	Reclining backrest (°)	-3+20	-3+20	-3+20	Eolo max load (kg)	25	35	50	Eolo weight (kg)	18	21	26
		Size 1	Size 2	Size 3																																										
<b>EOLO</b>	Seat width (cm)	17-24	22-29	27-34																																										
	Seat depth (cm)	19-26	22-34	31-43																																										
	Trunk height (cm)	26-38	33-45	41-53																																										
	Trunk width (cm)	bracket A 19-25	bracket A 22-28	bracket A 27-33																																										
		bracket B 14-20	bracket B 17-23	bracket B 22-28																																										
	Distance Popliteal Cable-Footboard (cm)*	10-21	10-26	10-36																																										
	Armrest height (cm)	14-19	20-26	20-26																																										
	Reclining backrest (°)	-3+20	-3+20	-3+20																																										
	Eolo max load (kg)	25	35	50																																										
Eolo weight (kg)	18	21	26																																											

<p><b><u>Technical features</u></b></p>	<p><b>MATERIALS:</b></p> <ul style="list-style-type: none"> <li>- metal structure made of aluminum and/or iron</li> <li>- the padding made of polyurethane foam of different density, bearing capacity and shape memory;</li> <li>- covers are made of PU/PES combined materials.</li> </ul> <p><b>MANTEINANCE AND CLEANING:</b> The removable cover can be washed at 30° with neutral detergent and sanitized with the following detergents/disinfectants:</p> <ul style="list-style-type: none"> <li>- PERSIL EXPERT coldzyme/HENKEL</li> <li>- 70% etanol</li> <li>- PERFORM/Schuelke</li> <li>- TPH PROTECT/Schuelke</li> <li>- MIKROZID AF liquid/Oktal pharma</li> <li>- DESCOGEN liquid/Antiseptica</li> <li>- MANORAPID/Antiseptica</li> <li>- PLIVASEPT/Pliva</li> <li>- INCIDIN PLUS/Ecolab</li> </ul> <p>Non usare detergenti a base di ipoclorito di sodio.</p> <p><b>PACKAGING:</b> Comes inside a carton</p>
<p><b><u>Accessories</u></b></p>	<p>Folding contouring thoracic supports  Dynamic butterfly harness  Pelvic belt  Contoured pelvic belt  Dynamic abduction harness  Confort headrest  Padded abductor wedge  Padding for footboard  Padding for tray  Head cover  Seat cover  Armrest cover (2pcs)  Thoracic support cover (2pcs)</p>
<p><b><u>Product image</u></b>  <i>(images are purely illustrative)</i></p>	