

## SENSO INPOCKET Positioning Backrest Technical Sheet

<b><u>Manufacturer identification</u></b>	<p>Pro Medicare s.r.l. Via Montagna Z.I. Lotto 41 72023 Mesagne (Br) P.IVA 01803920741 tel. 00 39 0831 777840 website: <a href="http://www.promedicare.eu">www.promedicare.eu</a> e-mail: info@promedicare.it</p>
<b><u>Device Commercial Name</u></b>	SENSO INPOCKET
<b><u>UDI-DI BASIC</u></b>	805660010POSTURAVERSAS3
<b><u>Device code</u></b>	<b>D05-701-X</b> (where X refers to the size; X=1-> refers to size 3630; X=2-> size 3635; X=3-> size 4030; X=4-> size 4035; X=5-> size 4430; X=6-> size 4435; X=7-> size 4830; X=8-> size 4835; X=9 -> size 3640; X=10 -> size 4040; X=11 -> size 4440; X=12 -> size 4840).
<b><u>Code of Classification (CND)</u></b>	<b>Y1809</b>
<b><u>Medical Device Class</u></b>	SENSO INPOCKET Positioning Backrest is CE marked according to EU Regulation 2017/745 and it is a non sterile Class I Medical Device.
<b><u>Intended use</u></b>	<p>The modular postural backrest for wheelchair users consisting of:</p> <ul style="list-style-type: none"> <li>- a flexible base equipped with fastening systems (hook-and-loop fasteners)</li> <li>- a stiffening system made up of adjustable-length tensioning straps equipped with fastening systems (hook-and-loop fasteners)</li> <li>- two padded and covered modular elements, equipped with fastening systems (hook-and-loop fasteners) with inside, a series of inserts with various shapes designed to accommodate and/or follow the shape of the user's rachis</li> </ul> <p>The SENSO INPOCKET backrest switches from an unfolded use configuration to a folded non-use configuration without having to be removed from the wheelchair. It is made from materials appropriately validated and tested to ensure lightness, comfort, and long-term durability, and it is also designed and manufactured in compliance with safety standards as required by the relevant Regulation (EU) 2017/745.</p> <p>Its intended use is to:</p> <ul style="list-style-type: none"> <li>- customize spinal posture;</li> <li>- manage postural disorders;</li> <li>- promote adequate positioning and support of the user's rachis;</li> <li>- promote distribution of body pressure;</li> <li>- relieve pressure on and ventilate the spine</li> </ul>
<b><u>Technical features</u></b>	<p><b>STANDARD CONFIGURATION:</b> the backrest comes complete with a flexible base element, a stiffening system made up of adjustable tensioning straps, two padded and covered modular elements with inside a series of inserts with various shapes designed to accommodate and/or follow the shape of the user's rachis.</p> <ul style="list-style-type: none"> <li>- <b>DIMENSIONS:</b> available in four sizes, each available in three different heights:</li> <li>- <b>3630-3635-3640</b> backrest with a width ranging from 35 cm to 38 cm (external tubes width), available in heights of 30 cm, 35 cm, and 40 cm.</li> <li>- <b>4030-4035-4040</b> backrest with a width ranging from 39 cm to 42 cm (external tubes width), available in heights of 30 cm, 35 cm, and 40 cm.</li> <li>- <b>4430-4435-4440</b> backrest with a width ranging from 43 cm to 46 cm (external tubes width), available in heights of 30 cm, 35 cm, and 40 cm.</li> <li>- <b>4830-4835-4840</b> backrest with a width ranging from 47 cm to 50 cm (external tubes width), available in heights of 30 cm, 35 cm, and 40 cm.</li> </ul> <p>The above-listed sizes are available in the following versions:</p> <ul style="list-style-type: none"> <li>- with 5 cm single straps;</li> <li>- with a 10 cm central strap and 5 cm single straps.</li> </ul> <p><b>WEIGHT:</b> &lt; 800 gr</p>

**Technical features**

**PERFORMANCE AND DURABILITY:** The expected duration of safe and effective performance of the devices is 2 years (padding) valid under normal conditions of use and proper maintenance in accordance with the specifications of the instruction manual.

**MATERIALI:** It is made of high-performance materials widely used in the medical industry:

1. the paddings made of polyurethane foam of different density, bearing capacity and shape memory and is covered with a material whose characteristics are listed in the following section;
2. Covers and the flexible base are made of PU/PES combined materials with the following characteristics:
  - antibacterial
  - non citotoxic
  - non irritating
  - resistant to virus penetration
  - antifungal
  - fireproof
  - biocompatible
  - OEKO-TEX®
  - Latex free
  - breathable
  - waterproof
3. The inserts are made of soft-density polyethylene material
4. The stiffening system, consisting of adjustable tensioning straps, is made of polymeric material

**MANTEINANCE AND CLEANING:** The removable cover can be washed at high temperatures 60° and sanitized with the following detergents/disinfectants:

- PERSIL EXPERT coldzyme/HENKEL
- 70% etanol
- PERFORM/Schuelke
- TPH PROTECT/Schuelke
- MIKROZID AF liquid/Oktal pharma
- DESCOGEN liquid/Antiseptica
- MANORAPID/Antiseptica
- PLIVASEPT/Pliva
- INCIDIN PLUS/Ecolab

Do not use detergents containing sodium hypochlorite. Do not spin-dry. Do not wash the structural parts; the padding covers can be superficially sanitized with the products listed above.



**PACKAGING:** Supplied inside a catron box

**Product image**  
(images are purely illustrative)

