




▶ **Multi-Functional Tilting Highchair
with High Low Indoor Base**

*For home, school, rehabilitation centers,
communities.....*


**CLASS I MEDICAL DEVICE
EU REGULATION 2017/745**





Reclining and Height adjustable Backrest

Length Adjustable Seat



Flip, Multi-Adjustable and Elevating Footboard

Standard Configuration



Purple/Taupe

Colors

Blue/Pastel Green



Covers

Washable
Disinfectable
Waterproof
Breathable

Seat Tilt
-5°+30°

Some Accessories



Multi-Adjustable
Headrest



Butterfly Harness

Pelvic Belt



Technical Details



Eolo	Size 1	Size 2	Size 3
Seat Width (cm)	17-24	22-29	27-34
Seat Length (cm)	19-26	22-34	31-43
Trunk Height (cm)	26-38	33-45	41-53
Trunk Width (cm)	Hdw A 19-25	Hdw A 22-28	Hdw A 27-33
	Hdw B 14-20	Hdw B 17-23	Hdw B 22-28
Distance Popliteal Cable – Footboard (cm)	10-21	10-26	10-36
Armrest Height (cm)	14-19	20-26	20-26
Reclining Backrest (°)	-3+20	-3+20	-3+20
Max Load Eolo (kg)	25	35	50
Weight Eolo (kg)	18	21	26

edlo



MADE IN



info@promedicare.it