

ORDER FORM

Invoicing Details	Shipping Details
Company name : _____ Delivery date required : ___ / ___ / ___ Address : _____ Postal Code .: _____ Town : _____ Tel: _____ Fax: _____ Purchase order: _____ Date of the order: ___ / ___ / ___	Ship to : _____ _____ <p style="text-align: center;">CUSTOMER STAMP AND SIGNATURE</p>

Pro Medicare srl grants **24 months** warranty on production defects from the date of the first put into service and **12 months** on covers and components subjected to wearing. Warranty is void if the product is used in any other way than stated in the Instruction Manual.
 The images are purely indicative and may not reflect the standard configuration


Standard Configuration *Ermes-Tipco*



- Seat adjustable in width, depth ,asymmetry, adduction, abduction
- Kit of inserts for a customised alignment and positioning of the pelvis
- Reclining and height adjustable backrest
- Kit of inserts for a customised alignment and positioning of the trunk
- Multi-adjustable and removable containment and positioning kit:
 - pelvis (adjustable in width depth, height, adduction/abduction)
 - arm support (adjustable in width, depth, height, angle, intra/extra rotation)
- Male interface plate
- Tipco indoor frame adapted with female interface plate :
 - * Height adjustable seat by gas spring
 - * Tilt in space by gas spring from -5° up to ° +30°
 - * 4 swivel castors (2 with locks)
 - * Multi-adjustable push handle



Sizes of the user



1 RH	depth of the seat RH	_____
1 LH	depth of the seat LH	_____
2 RH	length of the leg RH	_____
2 LH	length of the leg LH	_____
3 RH	Depth of the thorax RH	_____
3 LH	Depth of the thorax LH	_____
4	Height of the shoulders from the seat plane	_____
5	Width of the thorax	_____
6	Width of the pelvis	_____

Ermes

- **2** colour ○ blue/grey ○ fuchsia/grey
- **3** colour ○ blue/grey ○ fuchsia/grey
- **4** colour wine/grey

Notes: ○ = Optional without up-charge □ = Optional with up-charge

Type of the Seat

- **MIDI Seat**
Lateral/Anterior control
(standard configuration)




- **DEEP Seat**
Total Control (PTS)
(optional without up-charge)



Greater closure of the seat-backrest angle (≈ 5°)


Thoracic Supports

Flat Pads (values expressed in cm.)



Size .	1	2	3	4
H	7	7	8	10
L	7,5	8,5	10	12

Contoured Pads (values expressed in cm)



Size.	1	2	3	4
H	6	7	8	10
L	9,5	11	13	16

Side	Size	Bracket
<input type="checkbox"/> RH		<input type="checkbox"/> Fixed <input type="checkbox"/> Swing-Away type*: <input type="checkbox"/> A <input type="checkbox"/> B
<input type="checkbox"/> LH		<input type="checkbox"/> Fixed <input type="checkbox"/> Swing-Away type*: <input type="checkbox"/> A <input type="checkbox"/> B


Side	Size	Bracket
<input type="checkbox"/> RH		<input type="checkbox"/> Fixed <input type="checkbox"/> Swing-Away type*: <input type="checkbox"/> A <input type="checkbox"/> B
<input type="checkbox"/> LH		<input type="checkbox"/> Fixed <input type="checkbox"/> Swing-Away type*: <input type="checkbox"/> A <input type="checkbox"/> B

* In the technical information at the end of this form you can find the dimensions reachable with the two types of brackets

Belts and Harnesses (values expressed in cm, see details in the technical information)

Chest pad
(with thoracic supports only)

Size	HxL
<input type="checkbox"/> 1	(8x10)
<input type="checkbox"/> 2	(7x13)
<input type="checkbox"/> 3	(10x13)



Abduction Harness

<input type="checkbox"/>	XS (2)*
<input type="checkbox"/>	S (3)*
<input type="checkbox"/>	M (4)*



Pelvic Belt

Size .	HxL
<input type="checkbox"/> 1	(6x11)
<input type="checkbox"/> 2	(7x13)
<input type="checkbox"/> 3	(7x15)



Dynamic Butterfly Harness

<input type="checkbox"/>	XS (2)*
<input type="checkbox"/>	S (3)*
<input type="checkbox"/>	M (4)*



Jacket Harness

<input type="checkbox"/>	XS (2)*
<input type="checkbox"/>	S (3)*
<input type="checkbox"/>	M (4)*



Contoured Pelvic Belt




Mis.	HxL
<input type="checkbox"/> 1	(6x11)
<input type="checkbox"/> 2	(8x15)

Value H is referred to the height of the central part of the pad





* the number in the brackets is referred to the relevant Ermes size. Therefore it is **advisable** to combine the default size.


Headrest with multi-axial joint and multi-adjustable arm

<input type="checkbox"/> Comfort <input type="checkbox"/> fixed hardware <input type="checkbox"/> detachable hardware		<input type="checkbox"/> Anatomic <input type="checkbox"/> fixed hardware <input type="checkbox"/> detachable hardware		<input type="checkbox"/> 3 Piece pads supportive Headrest <input type="checkbox"/> fixed hardware <input type="checkbox"/> detachable hardware	
	<input type="checkbox"/>	S		<input type="checkbox"/>	XS
	<input type="checkbox"/>	M		<input type="checkbox"/>	S
	<input type="checkbox"/>	G		<input type="checkbox"/>	M
					<input type="checkbox"/> 1 (7x13) <input type="checkbox"/> 2 (9x20)
(values expressed in cm of the occipital pad only)					


Other Accessories

	<input type="checkbox"/> Knee elevating footplate hanger <input type="checkbox"/> RH <input type="checkbox"/> LH €
	Footrest Options: Separate footplates: <input type="checkbox"/> RH <input type="checkbox"/> LH Slotted Separate footplates: <input type="checkbox"/> RH <input type="checkbox"/> LH <input type="checkbox"/> Footboard (it does not allow asymmetry, adduction and abduction)

	<input type="checkbox"/> Shoes Holders
---	---



	<input type="checkbox"/> Calf Supports <input type="checkbox"/> RH <input type="checkbox"/> LH
	<input type="checkbox"/> Calf Panel (it does not allow asymmetry, adduction and abduction)
	<input type="checkbox"/> Lateral Supporting Pad for legs <input type="checkbox"/> RH <input type="checkbox"/> LH
	<input type="checkbox"/> Calf Strap anterior/posterior padded

<input type="checkbox"/> Padding for Separate footplates <input type="checkbox"/> Padding for footboard
--

<input type="checkbox"/> Detachable Tray grey colour with body contouring recess and edges <input type="radio"/> standard version <input type="radio"/> larger version for daily activities	<input type="checkbox"/> Detachable Tray transparent with body contouring recess and edges <input type="checkbox"/> Padding <input type="checkbox"/> superior <input type="checkbox"/> inferior																
Note: The trays are matched to the relevant Ermes size and are not interchangeable																	
	<table border="1"> <thead> <tr> <th>cm</th> <th>L</th> <th>P</th> <th>P1</th> </tr> </thead> <tbody> <tr> <td>2</td> <td>40</td> <td>35</td> <td>12</td> </tr> <tr> <td>3</td> <td>48</td> <td>40</td> <td>13.5</td> </tr> <tr> <td>4</td> <td>58</td> <td>45</td> <td>15.5</td> </tr> </tbody> </table>	cm	L	P	P1	2	40	35	12	3	48	40	13.5	4	58	45	15.5
cm	L	P	P1														
2	40	35	12														
3	48	40	13.5														
4	58	45	15.5														

<input type="checkbox"/> Elbow lateral-posterior support <input type="radio"/> RH <input type="radio"/> LH

Arm Supports and hip guides **Sizes** . **Notes** : Such components match the relevant Ermes size and are not interchangeable

Contoured Hip Guide					Flat Arm Support				
	cm	2	3	4		cm	2	3	4
	H	6.5	7	8		H	7.5	7.5	8.5
	L	17	20	23		L	17	20	23

Additional Outdoor Frame Adacta Giò for Ermes

Standard Configuration

- Adacta Gio' titanium satin anodised frame which plates have the same colour of the cover
- Adjustable single handlebar with folding backrest tubes
- Reclining mechanism by predetermined plate 3°,10°,18°,
- Adjustment of the seat depth
- Continuous Tilt in space by drive lever: from 2° up to 35° (from 6° to 39°)¹
- ¹tilt in space can be up-graded of 4° with an increment of the starting inclination of the frame by operating on the positioning holes of the casters forks
- Fixed anti-tip system
- Parking Brake by lever
- Adjustment of the seat depth
- Height adjustment: High (49 cm from the floor) / low (43 cm from the floor)*
- * (with 300mm and 500mm wheels only)
- ø 150 mm poly grey castors
- ø 300 mm poly grey posterior wheels with fast extraction
- Solid grey base assembled with female interface plate



<input type="checkbox"/>			
Size	<input type="checkbox"/> 2 (combined with Ermes size 2)	<input type="checkbox"/> 3 (combined with Ermes size 3)	<input type="checkbox"/> 6 (combined with Ermes size. 4)
Back Posts	<input type="checkbox"/> Fixed Push Handles Kit <input type="checkbox"/> Telescopic Push Handles Kit with Lever Block (max adjustment 10 cm) <input type="checkbox"/> Reclining Adjustable and separated handlebar kit <input type="checkbox"/> Quick-release pins for easier back reclining		
Castors	<input type="checkbox"/> 150 mm (6") poly <input type="checkbox"/> 175 mm (7") poly <input type="checkbox"/> grey <input type="checkbox"/> black		
Posterior Wheels	<input type="checkbox"/> 300 mm (12") poly <input type="checkbox"/> grey <input type="checkbox"/> black		<input type="checkbox"/> 300 mm (12") pneumatic
	<input type="checkbox"/> 300 mm (12") poly with drum brake <input type="checkbox"/> grey <input type="checkbox"/> black		<input type="checkbox"/> 300 mm (12") poly with drum brake + lever <input type="checkbox"/> <input type="checkbox"/> grey <input type="checkbox"/> black
	<input type="checkbox"/> 300 mm (12") pneumatic with drum brake		<input type="checkbox"/> 300 mm (12") pneumatic with drum brake + lever
	<input type="checkbox"/> 400 mm (16") poly black (only 45 cm height)		<input type="checkbox"/> 400 mm (16") poly black with drum brake (only 45 cm height)
	<input type="checkbox"/> 500 mm (20")poly without handrim		<input type="checkbox"/> 500 mm (20") pneumatic without handrim
<input type="checkbox"/> handrims for 500 mm wheels	<input type="checkbox"/> 500 mm (20") pneumatic with drum brake		<input type="checkbox"/> 500 mm (20") pneumatic with drum + lever

Technical Information

		Size.2	Size.3	Size.4
ERMES	Width of the Seat (cm)	≤25	18-28	22-35
	Depth of the Seat (cm)	19-26	22-34	31-43
	Height of the trunk (cm)	26-38	33-45	41-53
	Width of the thorax (cm)	staffa A 19-25	staffa A 22-28	staffa A 27-33
		staffa B 14-20	staffa B 17-23	staffa B 22-28
	Max abduction of the leg (°)	16	16	16
	Max adduction of the leg (°)	16	16	16
	Popliteal Space to footrest Length (cm)* * if the elevating knee angle is provided	8-19	8-29	8-39
	Height of the arm support (cm)	14-19	20-26	20-26
	Backrest Reclining (°)	-3+20	-3+20	-3+20
	Load of the Ermes (kg)	25	35	50
Weight of the Ermes (kg)	5.2	8.1	12.1	

		Size.2	Size.2	Size.3
TIPCO	Seat to floor height cm)	63	63	63
	Width overall dimension	67	67	67
	Depth overall dimensions (cm)	80	80	88
	Tilt in space (°)	-5+30	-5+30	-5+30
	Tipco Load(kg)	50	50	50
	Tipco Weight (kg)	13.6	13.6	14
	Total Weight (Tipco+Ermes) (kg)	18.8	21.7	26.1