

Ermes-Tipco Order form



Invoicing Details	Shipping Details			
Company name:	Ship to:			
Delivery date required:/				
Delivery address:	CUSTOMER STAMP AND SIGNATURE			
P. Code.: Town:				
Tel.: Fax:				
# Order: Date of the Order://				

Pro Medicare srl grants 24 months warranty on production defects from the date of the first put into service and 12 months on covers and components subjected to wearing. Warranty is void if the product is used in any other way than stated in the Instruction Manual.

Note: Any request for variation in dimensions, materials and type compared to the standard does not have CE marking, therefore the professional user becomes the manufacturer and has the obligation to draw up the technical documentation of the device.

Images are purely illustrative and may not fully reflect reality.

Ermes-Tipco Standard Configuration

- Seat adjustable in width, depth also asymmetrically in adduction and abduction
- Kit of inserts for a customised positioning of the pelvis
- Reclining and height adjustable backrest
- Kit of inserts for a customised positioning of the trunk
- Multi-adjustable and removable containment and positioning kit:
 - pelvis (adjustable in width, depth, height, adduction/abduction)
 - arm support (adjustable in width, depth, height, angle, intra/extra rotation)
- Male interface plate
- Tipco indoor frame adapted with female interface plate:
 - Height adjustable seat by gas spring
 - Tilt in space by gas spring from -5° up to ° +30°
 - 4 swivel castors (2 with locks)
 - Multi-adjustable push handle

Sizes of the user



1 RH	depth of the seat RH	
1 LH	depth of the seat LH	
2 RH	length of the leg RH	
2 LH	length of the leg LH	
3 RH	depth of the thorax RH	
3 LH	depth of the thorax LH	
4	height of the shoulders from the seat plane	
5	width of the thorax	
6	width of the pelvis	

Ű

<u>Ermes</u>								
0	size 2	color	0	blue/grey	0	fuchsia/grey		
0	size 3	color	0	blue/grey	0	fuchsia/grey		
0	size 4	color	0	wine/grey				
Note: If Ermes' sizes 2/3 are choosen, Tipco indoor base will be in white color; if Ermes' size 4 is choosen, Tipco indoor base will be in grey color								
Note: O = optional without up-charge 🗆 = optional with up-charge								

Type of the seat

0

MIDI Seat

lateral/anterior control (standard configuration)



DEEP Seat total control (PTS) optional without up-charge

0



□ Greater closure of the seat-backrest angle (±5°)

Thoracic Supports □ Flat pads Contoured pads Size (cm) 1 2 3 4 Size (cm) 1 2 3 4 н 7 7 н 6 7 8 10 8 10 7.5 8.5 L 10 12 L 9.5 11 13 16 s= 2.5 cm s= 2.5 cm Bracket type* Side Size Bracket type* Bracket shape Side Size Bracket shape □ Fixed $\Box A$ □ Fixed $\Box A$ \Box RH □ RH □ Swing-Away □ Swing-Away □ Fixed □ Fixed □ Swing-Away □ Swing-Away *In the technical information at the end of this form you can find the dimensions reachable with the two types of brackets (A type = L/90° Bracket; B type = 2.5 cm Offset Bracket) Belts and Harnesses (values expressed in cm, see details in the technical information) Chest pad □ Dynamic Adbuction Harness □ Pelvic belt XS (2)* HxW HxW Size Size □ 1 (8x10) П S (3)* □ 1 (6x11) □ 2 (9x13) M (4)* □ 2 (7x13) (with thoracic supports (10x13) □ 3 (7x15) □ 3 only) Dynamic Butterfly Harness Contoured Pelvic Belt Jacket Harness XS (2)* XS (2)* Size HxW S (3)* S(3)* □ **1** (6x11) П П M (4)* M (4)* \Box □ 2 (8x15) Value H is referred to the height of the central part of the pad *The number in the brackets is referred to the relevant Ermes size. Therefore it is advisable to combine the default size. Hardrest with multi-axial joint and multi-adjustable arm Confort with Anatomic with 3 piece pads supportive headrest with □ detachable hardware □ detachable hardware detachable hardware □ fixed hardware □ fixed hardware □ fixed hardware S XS (7x13) П □ 1 \Box М S □ 2 (9x20) values expressed in cm of G Μ П П the occipital pad only Other accessories Knee elevating legrest \Box RH 🗆 LH (with cover for angle knee)

separate footplates:

slotted separate footplates:

Footplate options:

□ footboard (it does not allow asymmetry, adduction and abduction)

		Shoes Holders							
Calf support Calf Support									
		□ Calf panel (it does no	ot allow asy	vmmetry,	adduction and	d abduction)			
		ateral supporting page	-	□ RH		· · · · · · · · ,			
R		Padded anterior/pd	-						
		Padding for separat		-					
a series		 Padding for footboa 	-						
Elbow lateral-pos		-							
Grey detachat	olo Trav				stan	ndard version dimens	sions		
-	ouring recess and	edges			6	Size (cm)	L	Р	P1
O standard ver	-	cubco			c.	2 (for Ermes size 2)	40	35	12
	n for daily activit	ies	-	Č.			40	40	12 5
	,,				P	3 (for Ermes size 3)	48	40	13.5
🗆 Transparent d	etachable Tray	,		P	J . I	4 (for Ermes size 4)	58	45	15.5
with body conto	ouring recess and	edges			larger version	on for daily activities	dimens	ions	
(only standard ver	sion)				6	Size (cm)	L	Р	P1
					c.	1 (for Ermes size 2)	50	45	13,7
		<u>tandard</u> grey color version ne, the naming changes (th		C.		2 (for Ermes size 3)	60	50	15.4
		ransparent tray, 3->2, 4->3			Pi				
				P V	<u> </u>	3 (for Ermes size 4)	68	55.6	19.3
Padding of the	tray 🗆 upp	er 🗆 lower							
Addictional	Outdoor Fra	me Adacta Giò fo	r Frmes						
 Standard Configuration Adacta titanium satin anodised frame which plates have the same color of the cover Adjustable single handlebar with folding backrest tubes Reclining mechanism by predetermined plate w/pin 3°, 10°, 18° Continuous Tilt in space by drive lever: from 2° up to 35° (from 6° up to 39°) tilt in space can be up-graded of 4° with an increment of the starting inclination of the frame by operating on the positioning holes of the casters forks Fixed anti-tip system Wheels Kit with brake Height adjustment: high (49 cm from the floor)/low (43 cm from the floor) with 300mm and 500mm wheels only Solid grey base assembled with female interface plate 									
		out up-charge							
Frame size	Z (combined	with Ermes size 2)	3 (con		ith Ermes siz	·	ned wit	n Ermes s	size 4)
<u>Castors</u>	🛇 150mm (6")	poly grey	O black		O 175mm (7") poly O grey	,	O black	
Ruote 300mm (12") poly grey with drum brake and parking brake lever 0 300mm (12") poly black with drum brake 300mm (12") poly black with drum brake 300mm (12") poly grey with drum brake 300mm (12") poly black with drum brake 300mm (12") poly black with drum brake and parking brake lever 300mm (12") poly grey with drum brake and parking brake lever 300mm (12") poly grey with drum brake and parking brake lever 1 300mm (12") poly grey with drum brake and parking brake lever 1 300mm (12") poly grey with drum brake and parking brake lever 1 300mm (12") poly grey with drum brake and parking brake lever 1 400mm (16") poly black with parking brake lever (only 45 cm height)									
	□ 400mm (16") poly <i>black</i> with drum brake (only 45 cm height)								
□ 400mm (16") poly <i>black</i> with drum brake and parking brake lever (only 45 cm height)									
) poly grey without hand				5 11			
) pneumatic <i>grey</i> without hand		-		۶r			
		r 500mm (20") wheels		, with park	ING DIAKE IEVE	-1			
1			on defense suit	+h d	aka				
) pneumatic <i>grey</i> with h							
□ 500mm (20") pneumatic <i>grey</i> with handrims with drum brake and parking brake lever									

	O Kit for fixed push handles							
	Kit for telescopic and angle adjustable single handlebar							
Back Posts	\square Kit for telescopic push handles with lever block (max adjustment 10 cm)							
	Kit for angle adjustable separated handlebar							
	Quick-release pins for easier back reclining							
<u>Tilt in space</u>	vith pedal O with lever							
Canopy adaptor included 🗆 Blue 🗆 Pink 🗆 Bordeaux								
Basket								
Rain cover								
<u>Foot muff</u>	🗆 Blue 🗆 Pink 🗆 Bordeaux							

TECHNICAL INFORMATION

Ermes		Size 2	Size 3	Size 4
Width of the Seat (cm)	≤25	18-28	22-35	
Depth of the Seat (cm)	19-26	22-34	31-43	
Height of the trunk (cm)	26-38	33-45	41-53	
		bracket A 19-25	bracket A 22-28	bracket A 27-33
Width of the thorax (cm)	bracket B 14-20	bracket B 17-23	bracket B 22-28	
Max abduction of the leg (°)	16	16	16	
Max adduction of the leg (°)		16	16	16
Popliteal Space to footplate Length (cm)* *if the elevating knee angle is provided		8-19	8-29	8-39
Height of the arm support (cm)		14-19	20-26	20-26
Backrest Reclining (°)		-3+20	-3+20	-3+20
Ermes load (kg)		25	35	50
Ermes weight (kg)		5.2	8.1	12.1
Contoured hip guide		Size 2	Size 3	Size 4
				5120 4
	Н	6.5	7	8
p	Р	17	20	23
Flat arm support		Size 2	Size 3	Size 4
	L	7.5	7.5	8.5
P	Р	17	20	23
Note: the postural unit for armrest and the lateral hip pad are	e comb	ined to the relative size o	of Ermes and are not inte	rchangeable
Тірсо		Size 2	Size 2	Size 3
Seat to floor height (cm)		63	63	63
Width overall dimension (cm)		67	67	67
Depth overall dimensions (cm)		80	80	88
Tilt in space (°)		-5+30	-5+30	-5+30
Tipco Load (kg)		50	50	50
Tipco Weight (kg)	13.6	13.6	14	
Total Weight (<i>Tipco+Ermes</i>) (kg)		18.8	21.7	26.1
Adacta Giò		Size 2	Size 3	Size 6
Seat to floor height (cm)		43	43	43
Width overall dimension (400/500 wheels) (cm)		48	53	58
Min/max depth overall dimensions (cm)		73-100	75-102	81-109
Tilt in space (°)		2-35 (from 6 up to 39) tilt in space can be up-graded of 4° with an increment of the starting inclination of the frame by operating on the positioning holes of the casters forks		
Predetermined backposts reclining (°)	3-10-18			
Adacta Giò load (standard configuration) (kg)	30	50	50	
Adacta Giò weight (kg)		12.8	13.4	14.5
Total weight (Adacta Giò+Ermes) (kg)		18	21.5	26.6



PRO MEDICARE S.r.l.

Via Montagna Z.I. Lotto 41-72023 Mesagne (Br)-ITALY Tel: +39-0831-777840 Fax:+39-0831-730739 Website: <u>www.promedicare.eu</u> e-mail: sales@promedicare.it